

CANCELLOUS GRANULES IN THE APPLICATOR

Cancelous cubes

Origin: Human
Tissue: Cancelous bone
Processing: Allotec® purification procedure
Inactivation: Min. SAL10⁻⁶ for viruses and bacteria
Sterilisation: Gamma irradiation
Application: Bone void filler
Rehydration: Min. 10 minutes

The applicator of C⁺TBA is a special form of primary packaging for cancellous bone granules, which simplifies both the rehydration with a physiological medium as well as the application of the granules into the defect zone.



ORDER INFORMATION*

DESCRIPTION	ITEM NUMBER	VOLUME [cc]
Cancelous granules in the applicator	ALO360	7
	ALO361	15
	ALO362	30



* Please note: Due to the nature of human bone tissue and the technical possibilities of shaping, slight deviations of the specified sizes may occur.

BT Medicals AB
nordics & baltics

Scheelevägen 30-32
223 63 Lund

+46 364404444
office@btmedicals.com

www.btmedicals.com