


CANCELLOUS RING

The cancellous ring was specially developed to simplify the reconstruction of the glenoid in a total shoulder prosthesis.

cancellous ring

Origin:	Human
Tissue:	Cancellouse bone
Processing:	Allotec® purification procedure
Inactivation:	Min. SAL10 ⁻⁶ for viruses and bacteria
Sterilisation:	Gamma irradiaton
Application:	Remodelling of the glenoid in case of shoulder endoprosthesis
Rehydration:	Min. 10 minutes



ORDER INFORMATION*

DESCRIPTION	iD [mm]	ITEM NUMBER	D [mm]	H [mm]
Cancellous ring	1,5	ALO431	26	10
		ALO433	32	10
		ALO432	26	20
		ALO434	32	20
	7,7	ALO436	26	10
		ALO437	32	10
		ALO435	26	20
		ALO430	32	20



* Please note: Due to the nature of human bone tissue and the technical possibilities of shaping, slight deviations of the specified sizes may occur.

BT Medicals AB
nordics & baltics

Scheelevägen 30-32
223 63 Lund

+46 364404444
office@btmedicals.com

www.btmedicals.com