


# CANCELLOUS CYLINDER

Cancellous bone cylinders are preferably used in sports medicine for filling drill channels in cruciate ligament revisions.

## Cancellous cylinder

Origin: Human  
Tissue: Cancellous bone  
Processing: Allotec® purification procedure  
Inactivation: Min. SAL10<sup>-6</sup> for viruses and bacteria  
Sterilisation: Gamma irradiation  
Application: Tunnel filling  
Rehydration: Min. 10 minutes



## ORDER INFORMATION\*

DESCRIPTION	ITEM NUMBER	D [mm]	H [mm]
Cancellous cylinder	ALO423	10	20
	ALO424	10	30
	ALO425	12	20
	ALO426	12	30
	ALO427	14	20
	ALO428	14	30



\* Please note: Due to the nature of human bone tissue and the technical possibilities of shaping, slight deviations of the specified sizes may occur.

**BT Medicals AB**  
nordics & baltics

Scheelevägen 30-32  
223 63 Lund

+46 364404444  
office@btmedicals.com

www.btmedicals.com