

AlloHeal® C+TBA

DBM PUTTY

AlloHeal® C+TBA DBM Putty consists of 100% demineralised, allogenic bone of human origin with no carrier substances. AlloHeal® C+TBA DBM Putty is proven to be osteoinductive, meaning it promotes bone growth, and thereby contributes to rapid bone healing.

AlloHeal® C+TBA DBM Putty

Source: Human
Tissue: 100% allogenic bone tissue
Preparation: Allotec® purification procedure & demineralisation
Inactivation: Min. SAL 10⁻⁶ for viruses and bacteria
Sterilisation: Gamma irradiation
Preservation: Ready to use
Application: Repairing bone defects



ORDER INFORMATION

DESCRIPTION	ITEM NUMBER	GTIN	V [cc]
AlloHeal® C+TBA DBM Putty	81010	9009808006505	1,0
	81025	9009808006512	2,5
	81050	9009808006529	5,0
	81100	9009808006536	10,0



BT Medicals AB
nordics & baltics

Scheelevägen 30
223 63 Lund

+46 364404444
office@btmedicals.com

www.btmedicals.com